

ALL-CERAMIC DENTAL IMPLANT SOLUTIONS



CeraRoot

ZIRCONIUM OXIDE DENTAL IMPLANTS


CE

FDA

 Health Canada  Santé Canada

Scientific Evidence

Bone-to-implant contact of 78% at 3 months.
One piece implant = no prosthetic connections



References

1. One-year follow-up of first consecutive 100 zirconia dental implants in humans: a comparison of 2 different rough surfaces. *Int J Oral Maxillofac Implants*. 2007 May-Jun;22(3):430-5.
2. Erste klinische Einjahresergebnisse von 100 Zirkonoxidimplantaten. Ein Vergleich zwischen zwei unterschiedlich rauhen Oberflächentypen. *Implantologie* 2007;15(4):429-436.
3. Zirconia implants and all-ceramic restorations for the aesthetic replacement of two central incisors. *EJED* 2008;3(2):175-185.
4. Ovoid Zirconia Implants: Anatomic Design for Premolar Replacement. *Int J Period Rest Dent* 2008; 28:609-615.
5. Ovale Zirkonoxidimplantate: Anatomisches Design als Ersatz für Prämolare. *Int J Par Rest Zahn* 2008;28:583-589.
6. Replacement of congenitally missing maxillary permanent canine with a zirconium oxide dental implant and crown. A case report from an ongoing clinical study. *Oral Surgery* 2008;1(2):140-144.
7. Full mouth oral rehabilitation in a titanium allergy patient using zirconium oxide dental implants and zirconium oxide restorations A case report from an on-going clinical study. *EJED*; 2010 Summer; 5(2):190-203.
8. Five year success rate of first consecutive 831 zirconia dental implants in humans. A comparison of three different rough surfaces. *JOMI* 2010 Mar-Apr;25(2):336-44.

Full articles available at: www.ceraroot.com/store/download/scientific-publications/

Company

Oral Iceberg is a global leader in metal-free, holistic and esthetic implant dentistry. Having pioneered the latest research and developments in the field of zirconium oxide implant dentistry with the CeraRoot system, we have a philosophy of doing more to advance dental regeneration, restoration, and replacement, as well as patient care.

With the CeraRoot implant system and CeraCrown prosthetic components, we deliver superior solutions that allow dental professionals to provide the best possible care to patients.

Located in Barcelona, Spain, ORAL ICEBERG company is a global leader in zirconium oxide implant and restorative dentistry. Founded by Drs. Oliva in 1995, the company has experienced a great growth with the development of the CeraRoot zirconium oxide implant system and the CeraCrown prosthetic components.

In collaboration with The Leading Dental Centers, research institutes and universities, we study and develop implants, prosthetic components and instruments for use in tooth replacement solutions.

We also offer comprehensive training and services to the dental professionals worldwide, including training and education, which is provided in collaboration with the LDCW group. Our products are available worldwide through the online shop.

Oral Iceberg encourages healthy living and social responsibility.

*The Leading Dental Centers
of The World®*



oral iceberg

ORAL ICEBERG, S.L.
Josep Umbert, 128 - 08402 Granollers
Barcelona - SPAIN
Phone +34 938 614 841
Fax +34 938 792 373
www.ceraroot.com
info@ceraroot.com
Skype id: oral.iceberg

CE₁₂₇₅

ISO 13485

Our products are conformance with European Union regulations.

FDA K093595

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician / dentist

 Health Canada Santé Canada **87566**


Health Canada



Implants


CeraRoot 12

FOR LATERAL OR LOWER INCISORS




CeraRoot 21

FOR SMALL CENTRAL INCISORS AND CANINES




CeraRoot 14

FOR PREMOLARS AND BICUSPIDS





CeraRoot 11

FOR CENTRAL INCISORS AND CANINES



CeraRoot 16

FOR MOLARS



Indications



CeraRoot 12



CeraRoot 21





CeraRoot 11



CeraRoot 14



CeraRoot 16





Osseointegration



The osseointegration that results from the enhanced bone formation and firm anchorage on Ice Surface is characterized by a remarkable bone-to-implant contact. Early loading is possible on CeraRoot implants when excellent stability is achieved in healthy patients with good bone quality.

The unique properties of the Ice Surface accelerate the osseointegration process and achieve secondary stability sooner than machined or blasted surface ceramic implants. This means that early loading protocols can be achieved with more predictability thanks to earlier osseointegration. The latest clinical results on Ice Surface confirm exceptional maintenance of marginal bone level after 5 years.

Ice Surface is currently being evaluated in several centers involving more than 800 CeraRoot implants placed in 378 patients. For the patient, this means the benefit of more bone more rapidly and the reliability of a long-term treatment outcome. CeraRoot is the first generation acid etched ceramic implant surface with more than 5 years' clinical follow-up. The absence of a prosthetic connection due to its one-piece structure accelerates the healing process and preserves the surrounding bone.








Ice Surface™

Ice Surface is the first acid etched surface applied on a zirconium oxide ceramic implant.

This new surface helps achieve optimal conditions for the integration between implant and bone and soft tissues. It is a rough surface especially designed for achieving faster and stronger osseointegration compared to machined zirconium oxide ceramic implants. The enhanced osseointegration and firmer anchorage in the surrounding tissues achieved with Ice Surface™ increase the predictability of implant treatment, especially in clinically demanding situations.



Mechanical Testing

The strength and machining characteristics of yttrium reinforced - zirconium oxide (γ -tzp) implants make this material ideal for the fabrication of durable, high-precision and highly esthetic ceramic implants and crowns.

The material properties, in combination with the design and production process technology of the CeraRoot implant system, result in an excellent mechanical stability.

To test the mechanical properties of the CeraRoot implant, fatigue strength testing according to the ISO 14801 standard was conducted.

With 5 million validation cycles at 10 Hz at the same benchmark load as the titanium implants, the CeraRoot implant met all acceptance criteria. The angle of loading force was at 30° to the axis of the implant body (Fig.1).

This test result show that the mechanical properties offered by implant type CeraRoot 11 indicated for Central Incisors and Canines has an excellent mechanical properties even in the most severe conditions.

The simulation of 5.000.000 cycles is to simulate an oral environment function of about 20 years (Fig. 2).

Zirconia has a flexural strength and fracture toughness almost twice as high as that of alumina, which makes zirconia very resistant to masticatory forces, while still maintaining an exact precision.

These proven mechanical qualities allow the restorative team high flexibility for indication and application.

Partially-stabilized zirconia powder with uniform dispersion of 3 mol % yttria exhibits superior sintering properties and higher aging resistance at a lower sintering temperature of 1350 °C.

Sintered bodies produced with primary powder show a fine crystal grain structure resulting in great improvements in strength, fracture toughness, as well as resistance to wear and aging.

Zirconium Oxide has numerous applications such as materials for medical devices, industrial parts and everyday products.

CeraRoot implants are manufactured with compositions according to ISO 13356.

FIG. 1 - 30° Off-axis Fatigue Testing




FIG. 2 - CeraRoot 11, Fatigue Test results

| Test with CeraRoot 11 | Load (N) | N°. of cycles |
|-----------------------|----------|---------------|
| Static | | |
| Mean of 3 implants | 2730 | 1 |
| Dynamic | | |
| Mean of 8 implants | 1967 | 5.000.000 |



CeraCrown™

the crown is available before the implant is inserted

All CeraRoot implants have been designed to avoid any preparation of the abutment and prosthetic shoulder part. That is why each implant has its own indication.

This makes zirconia stronger over time and provides a perfect fit between the implant and the crown.



The professional can grind or polish the crown to its final shape and send it to the lab to be finished with porcelain layers or over-pressing technique.

The CeraCrown can be used as:

Impression coping

Temporary crown

Zirconia crown core



CeraCrown 12

CeraCrown 21

CeraCrown 11

CeraCrown 14

CeraCrown 16

These products are intentionally designed to be easy to use in your daily practice. Combining the finest dental prosthetics with the world's only industrialized process customized crowns, we ensure you really can deliver functionality and esthetics to all your patients.

CeraRoot's unique CeraCrown™ technology offers esthetic and functional dental restorations for all teeth. Based on the latest manufacturing technologies, the CeraCrown system provides completely individualized prosthetics with unbeatably precise fit for crowns.

A combination of biocompatibility, beauty, and strength is guaranteed. Zirconia refracts and transmits light in much the same way as a natural tooth, thereby giving the restoration a natural looking appearance. In fact, the end result is often an improvement on nature. In thin periodontium type, the white body of the implant is never visible through the gum.

Beautiful gingiva. Best biomaterials, surface, design and procedures created for maintaining and regaining natural soft tissue. Supra or paragingival cementation of the crowns has never been easier.



Assembling of CeraRoot and CeraCrown



CeraRoot 16



CeraCrown 16




Crown finished

Drill Sequence

CeraRoot 12

FOR LATERAL OR LOWER INCISORS




- 1) pilot drill
 - 2) 2.3 mm
 - 3) 2.8 mm
 - 4) 3.5 mm
 - 5) 4.2 mm (only the first 3 mm)
- NO countersink!

ATTENTION: Do NOT exceed 35N of torque during insertion of implant

CeraRoot 21


FOR SMALL CENTRAL INCISORS AND CANINES



- 1) pilot drill
- 2) 2.3 mm
- 3) 2.8 mm
- 4) 3.5 mm
- 5) Countersink 21 ref. 004.021

CeraRoot 14


FOR PREMOLARS AND BICUSPIDS



- 1) pilot drill
- 2) 2.3 mm
- 3) 2.8 mm
- 4) 3.5 mm
- 5) Countersink 14 ref. 004.014

CeraRoot 11


FOR CENTRAL INCISORS AND CANINES



- 1) pilot drill
- 2) 2.3 mm
- 3) 2.8 mm
- 4) 3.5 mm
- 5) 4.2 mm
- 6) Countersink 11 ref. 004.011

CeraRoot 16

FOR MOLARS



- 1) pilot drill
- 2) 2.3 mm
- 3) 2.8 mm
- 4) 3.5 mm
- 5) 4.2 mm
- 6) Countersink 16 ref. 004.016

Instruments

Implants

| SURGICAL PARTS | |
|----------------|--|
| Reference | Description |
| 001.001 | Initial bur 2mm (triangle shaped) |
| 001.002 | Twist drill 2.3mm - steel |
| 001.003 | Twist drill 2.8mm - steel |
| 001.004 | Twist drill 3.5mm - steel |
| 001.005 | Twist drill 4.2mm - steel |
| 001.006 | Bur extension |
| 001.007 | Right-angle/Handpiece adapter 4x4 |
| 002.002 | Twist drill 2.3mm - ceramic |
| 002.003 | Twist drill 2.8mm - ceramic |
| 002.004 | Twist drill 3.5mm - ceramic |
| 002.005 | Twist drill 4.2mm - ceramic |
| 003.016 | Countersink 16 - steel |
| 003.011 | Countersink 11 - steel |
| 003.021 | Countersink 21 - steel |
| 003.014 | Countersink 14 - steel |
| 004.016 | Countersink 16 - ceramic |
| 004.011 | Countersink 11 - ceramic |
| 004.021 | Countersink 21 - ceramic |
| 004.014 | Countersink 14 - ceramic |
| 005.114 | Mallet for implant 14 |
| 005.214 | Impactor for implant 14 |
| 005.016 | Transporter 16 Long- steel |
| 005.017 | Transporter 16 Short- steel |
| 005.011 | Transporter 11- steel |
| 005.021 | Transporter 21- steel |
| 005.012 | Transporter 12- steel |
| 006.016 | Transporter 16 Long- zirconia |
| 006.017 | Transporter 16 Short- zirconia |
| 006.011 | Transporter 11- zirconia |
| 006.021 | Transporter 21- zirconia |
| 006.012 | Transporter 12- zirconia |
| 020.001 | Big Surgical cassette empty |
| 020.003 | Small Surgical cassette empty |
| 020.004 | Ceramic tray 50x75 for disposal of drills |
| 020.007 | Surgical Kit with steel drills/instruments |
| 020.008 | Surgical Kit with ceramic drills/instruments |
| 005.500 | Steel Tweezers with zirconia tips |
| 030.001 | Torque Wrench 35N - titanium |
| 030.002 | Adapter 4x4 SHORT for Torque wrench |
| 030.003 | Adapter 4x4 LONG for Torque wrench |
| 030.004 | Torque Wrench Kit (030.001+ 030.002 + 030.003) |
| 030.005 | NobelBiocare® adapter |
| 030.006 | Straumann® adapter |

| PROSTHETIC PARTS | |
|------------------|--|
| Reference | Description |
| 007.016 | CeraCrown 16 - zirconia coping white color |
| 017.016 | CeraCrown 16 - zirconia coping A1 color |
| 027.016 | CeraCrown 16 - zirconia coping A2 color |
| 007.011 | CeraCrown 11 - zirconia coping white color |
| 017.011 | CeraCrown 11 - zirconia coping A1 color |
| 027.011 | CeraCrown 11 - zirconia coping A2 color |
| 007.021 | CeraCrown 21 - zirconia coping white color |
| 017.021 | CeraCrown 21 - zirconia coping A1 color |
| 027.021 | CeraCrown 21 - zirconia coping A2 color |
| 007.014 | CeraCrown 14 - zirconia coping white color |
| 017.014 | CeraCrown 14 - zirconia coping A1 color |
| 027.014 | CeraCrown 14 - zirconia coping A2 color |
| 007.012 | CeraCrown 12 - zirconia coping white color |
| 017.012 | CeraCrown 12 - zirconia coping A1 color |
| 027.012 | CeraCrown 12 - zirconia coping A2 color |
| 008.116 | Lab Analogue 16 for CeraRoot 16 - zirconia |
| 008.111 | Lab Analogue 11 for CeraRoot 11 - zirconia |
| 008.121 | Lab Analogue 21 for CeraRoot 21 - zirconia |
| 008.114 | Lab Analogue 14 for CeraRoot 14 - zirconia |
| 008.112 | Lab Analogue 12 for CeraRoot 12 - zirconia |

| Reference | Description | | |
|--------------------|---|-----------|-------------------|
| CeraRoot 16 | FOR MOLARS | *a | c s |
| 160.010 | CeraRoot 16 L= 10 mm | 4.8 | - 6.5 - 8 |
| 160.012 | CeraRoot 16 L= 12 mm | 4.8 | - 6.5 - 8 |
| 160.014 | CeraRoot 16 L= 14 mm | 4.8 | - 6.5 - 8 |
| CeraRoot 11 | FOR UPPER CENTRAL INCISORS AND UPPER CANINES | | |
| 110.010 | CeraRoot 11 L= 10 mm | 4.8 | - 6 - 6.5 |
| 110.012 | CeraRoot 11 L= 12 mm | 4.8 | - 6 - 6.5 |
| 110.014 | CeraRoot 11 L= 14 mm | 4.8 | - 6 - 6.5 |
| CeraRoot 21 | FOR UPPER CENTRAL INCISORS AND CANINES | | |
| 210.010 | CeraRoot 21 L= 10 mm | 4.1 | - 4.8 - 6 |
| 210.012 | CeraRoot 21 L= 12 mm | 4.1 | - 4.8 - 6 |
| 210.014 | CeraRoot 21 L= 14 mm | 4.1 | - 4.8 - 6 |
| CeraRoot 14 | FOR BICUSPIDS/PREMOLARS AND LOWER CUSPIDS | | |
| 140.008 | CeraRoot 14 L= 8 mm | 3.5 | - 4.8 - 5 |
| 140.010 | CeraRoot 14 L= 10 mm | 3.5 | - 4.8 - 5 |
| 140.012 | CeraRoot 14 L= 12 mm | 3.5 | - 4.8 - 5 |
| 140.014 | CeraRoot 14 L= 14 mm | 3.5 | - 4.8 - 5 |
| CeraRoot 12 | FOR UPPER LATERAL AND LOWER INCISORS | | |
| 120.010 | CeraRoot 12 L= 10 mm | 4.1 | - 4.1 - 4.1 |
| 120.012 | CeraRoot 12 L= 12 mm | 4.1 | - 4.1 - 4.1 |
| 120.014 | CeraRoot 12 L= 14 mm | 4.1 | - 4.1 - 4.1 |

*a --> diameter of implant at apical level
 c --> diameter of implant at bone crest level
 s --> diameter of implant at shoulder level

Drawings of parts are available on the website



CERAROOT TRAINING

Visit www.ceraroot.com
for available courses

CeraRoot

 oral iceberg

ORAL ICEBERG, S.L.
Josep Umbert, 128 - 08402 Granollers
Barcelona - SPAIN
Phone +34 938 614 841
Fax +34 938 792 373
www.ceraroot.com
info@ceraroot.com
Skype id: oral.iceberg

CE₁₂₇₅

ISO 13485

Our products are conformance with
European Union regulations.

FDA K093595

Caution: Federal Law (USA) restricts
this device to sale by or on the order of a
physician / dentist

 Health Canada Santé Canada 87566

Health Canada